



Easy Armboard™ With Trigger CE

With the trigger design, simply pull the radiolucent release mechanism for full articulation and release it to lock. A spring lever accessible from both sides securely locks the armboard to the rail and eliminates the uncertainty of gravity locks.

#O-ABSLT Easy Armboard With Trigger (US)

#O-ABSLTE Easy Armboard With Trigger (Eur)

#O-ABSLTUK Easy Armboard With Trigger (UK)

**Eliminate Broken
Armboard Teeth!**



What's Required		OR	
What's Optional			



Easy Armboard™ Lift-to-Position CE

Lift the end of the board and rotate to the desired position. The patient's arm weight locks the board in place when lowered. A spring lever securely locks the armboard to the rail eliminating the uncertainty of gravity locks.



#O-ABSL Easy Armboard™ Lift-to-Position (US)

#O-ABSLE Easy Armboard™ Lift-to-Position (Eur)

#O-ABSLUK Easy Armboard™ Lift-to-Position (UK)

#O-ABSLD Easy Armboard™ Lift-to-Position (Denyer)

What's Required		OR	
What's Optional			



Dura-Board® Armboard

Components are enclosed in a durable, molded plastic body to protect against denting and breakage. Easy, no-tilt mounting, removal and repositioning. Pull-trigger provides easy 180° lateral adjustment with 11 positions. Rail clamp latch prevents inadvertent detachment. Lightweight (3.8 lbs. / 1.7 kg. without pad) and radiolucent. Integral storage and carrying handle.



Fits US Rails Only

#A-30003 Dura-Board® Armboard

Choose a Pad		OR	
What's Optional			



Easy Extra-Wide Armboard

The Extra-Wide Armboard features a 11" x 24" (28 cm x 61 cm) surface. Attaches easily and securely to the table rail, using our spring lever design. Triggerless design allows you to simply lift the end of the board and rotate to the desired position. When lowered, the patient's arm weight causes the board to lock securely.

Fits US Rails Only

#O-ABSLW Easy Extra-Wide Armboard

What's Included		
What's Optional		

**All Armboards
are Radiolucent!**