

C-Flex™ Head Positioning System

CE

Just Squeeze and Position!

The modular system can be used in all types of spine surgery. The quick-connect coupler makes connecting to a Mayfield skull clamp a snap! Ability to “squeeze and position” supports readjustments intraoperatively. An extensive range of motion makes C-Flex ideal for use with extremely kyphotic patients, anterior artificial discs, and lateral procedures.

- Cervical flexion / extension
- Lumbar positioning
- Adjust vertically for kyphotic patients
- Lateral positioning
- Device weight = 20 lbs. (9 kg)
- 6" (15 cm) lateral adjustment range
- 28" (71 cm) height adjustment range
- 9" (23 cm) extension adjustment range
- +/- 15° lateral tilt
- Works with the Mayfield® Skull Clamp, and Mayfield Horseshoe
- Includes fold-away mirror
- Eyes free from pressure
- Easy airway access
- 500 lbs. (227 kg) patient weight limit
- Fits the Allen Flex Frame
- Fits the Jackson Spine Table

#A-70700 C-Flex Head Positioning System



Anterior Access



Kyphotic Patients











Prone



Clean Cape Disposables



Just Squeeze and Position

What's Included	What's Required	What's Optional
 C-Flex Base Unit #A-70701  Comfort Mask Disposables #A-70310 Includes Case of 5  C-Flex Clean Cape Disposables #A-70702 Includes Case of 10	 Starburst Adaptor #A-70710 1 Included  Prone Mask Head Module #A-70715 1 Included  Flat Plate Head Module #A-70720 1 Included	 Reverse Starburst Adaptor #A-70725
	 Comfort Mask Disposables #A-70310 Case of 5  C-Flex Clean Cape Disposables #A-70702 Case of 10	