



For Rectangular Blades

For Round Posts

Tri-Clamp®

Very versatile clamp, three independent adjustments for lithotomy, abduction and rail securement. Attach the clamp anywhere on the OR table rail between rail standoffs. Accepts ½" to ⅝" (1.3 cm to 1.6 cm) round and 1" x ¼" (2.5 cm x 0.6 cm) flat accessories.

#A-40016 Tri-Clamp (US)

#A-40016E Tri-Clamp (Eur)



For Rectangular Blades

For Round Posts

Universal Clamp

CE

Single clamp secures accessories with round and rectangular mounting posts from any manufacturer. Accessories can be mounted either vertically or horizontally (e.g. foot extension) with the same clamp. Attach the clamp anywhere on the OR table rail between rail standoffs. Accepts ⅝" (1.6 cm) round and 1" x ¼" (2.5 cm x 0.6 cm) flat accessories. Not recommended for use with stirrups.

#O-RCP2 Universal Clamp (US/Eur)



For Rectangular Blades

For Round Posts

Universal Accessory Clamp

Unique, basic clamp that fits both ½" to ⅝" (1.3 cm to 1.6 cm) table accessories with round posts as well as blade-type posts ¼" x 1" (.6 cm to 2.5 cm). Attach the clamp anywhere on the OR table rail between rail standoffs. Not recommended for use with stirrups.

#A-40015 Universal Accessory Clamp (US Only)



For Round Posts

Direct Clamp

A strong and light-weight replacement for "Clark Sockets." Increased clamping force stops rotation of legholders. Teeth are concealed for easier cleaning. Unique design prevents accidental sliding along rail when not in use. Attaches at OR table rail ends or rail notches. Will not hang up on rail-end gravity locks. Indicator marks facilitate faster, easier setup. Accepts up to ⅝" (1.6 cm) round post.

#O-CS Direct Clamp (US)

#O-CSU Direct Clamp (Eur/UK)

#O-DCD Direct Clamp (Denyer)

#O-CSSW Direct Clamp (Swiss)



For Round Posts

Simple Clamp


CE

Innovative design allows placement under/over drapes. Durable aluminum and stainless steel construction. Attach the clamp anywhere on the OR table rail between rail standoffs. Holds all ½" to ⅝" (1.3 cm to 1.6 cm) diameter accessories.

#A-40018 Simple Clamp (US/Eur)

#A-40018-UK Simple Clamp (UK)

#A-40018-DEN Simple Clamp (Denyer)



For Rectangular Blades

Rail Clamp

CE

Securely locks all accessories with flat mounting blades anywhere along the table side rail. Attaches at OR table rail ends or rail notches. Perfect for tabletop extensions and knee/foot board applications. Strong and light-weight replacements for original rail locks. Accessories can be mounted either vertically or horizontally. Reduces the "free play" associated with standard clamps. Accepts 1" x ¼" (2.5 cm x 0.6 cm) flat accessories.

#O-RC2 Rail Clamp (US)

#O-RC2XUK Rail Clamp (Eur/UK)

#O-RC2XD Rail Clamp (Denyer)

#O-RC29X32MM Rail Clamp (Jpn)

#O-RCS Rail Clamp (Swiss)