



No Lifting Required!


### Patient Transfer Board

#### Roller-Free Roller Board - Better For Patients, Easier On Nurses, and Inexpensive

- **No Rollers** - Less traumatic transfer for patients
- **Comfortable** - Soft but sturdy board is easy on patients
- **Lightens the Load** - No lifting required, even with your heaviest patients
- **Lightweight** - Very easy to use
- **No Special Training Requirements** - Uses same or similar transfer technique
- **Radiolucent** - You can image with the board under the patient
- **Fowler-Position Transfer** - Long board allows for transfer in seated position
- **Foldable** - Long board folds for easy storage and set-up
- **Straps** - Carry this lightweight device on your shoulder and keep your hands free



Makes Fowler-Position Transfer a Snap!



Ultra-Lightweight!  
Short Board is 3 lbs. (1.5 kg)  
Long Board is 6.3 lbs. (3 kg)

#A-83000 Long Patient Transfer Board  
71" x 20" (180 cm x 50 cm)

#A-83001 Short Patient Transfer Board  
39" x 16" (100 cm x 40 cm)


#A-83002 Disposable Covers (Case of 50)  
(sized for Long Boards)



#A-83001



#A-83000



### Warrior Evacuation Litter Pad

The Warrior Evacuation Litter Pad is a variation of the 2" Astrus™ Stretcher pad made by Allen Medical Systems. The Astrus stretcher pad is used for enhanced tissue pressure management in the hospital arena. The Astrus pad has been modified to Air Force Air Evacuation specs and renamed the "Warrior Evacuation Litter Pad" aka "WELP." The WELP is designed as an enhanced tissue pressure management litter surface for use by the D.O.D in wartime medical evacuation missions. There are two pressure relief vents built into the pad located one at each end at opposing sides.

Dimensions: 72" x 23" x 2"

Weight: 10.5 lbs. per pad

National Stock Number - 6530-01-548-0262

#H-223-48-223-A1 Warrior Evacuation Litter Pad