



#A-70002 with C-Flex



#A-70000 with C-Prone

## Allen Spine System



### Give Your OR Table Jackson Table Capabilities!

*A Fraction of the Price of a Standard Spine Table*


The Allen Spine System is an extension for your OR table that enables you to perform spine surgery with the positioning benefits of a Jackson table! In addition it offers intra-operative lumbar flex and a small storage footprint. The advanced accessories adjust to fit the individual patient's anatomy.

- Intra-operative Flex – use the power of the table to flex and reposition the patient
- Radiolucent Extension 48" x 17" (1220 mm x 430 mm) – provides excellent C-Arm access
- Saves Space – stores in less than 3 sq. ft., (1 sq. m) and accessories can be stored on the frame
- Optimized pressure management - chest, hip, and head supports
- Versatile – prone, lateral, and supine positioning
- Range of Motion – provides 20° lateral tilt, trend and flex powered by the OR table
- Height – 25" to 38" (635 mm - 965 mm)
- Fits OR tables with outside rail width 20" - 24" (510 mm - 610 mm)
- 500 lbs. (227 kg) patient weight capacity
- XLIF kit available, see page 20

#A-70002 Spine System (Premium System with C-Flex)

#A-70000 Spine System (Standard System with C-Prone)

### Allen Spine System Components



- A.** C-Prone Head Positioner
- B.** Standard Wingset
- C.** Flex Frame
- D.** Flex Frame Armboard Rail
- E.** Flex Frame OR Table Clamp Set
- F.** Prone Arm Supports
- G.** Easy Lock Socket
- H.** Supine Top Pads
- I.** Supine Top Set

### Spine System #A-70002

### Spine System #A-70000

	Choose Head Positioner	#A-70002 Spine System	#A-70000 Spine System
What's Included (A-70000 + A-70002)	<p>C-Flex Head Positioner System #A-70700</p>	<p>C-Prone Head Positioner #A-70300</p>	
	<p>Flex Frame #A-70100</p>	<p>Prone Arm Supports #A-70610 Includes 2</p>	<p>Supine Top Set #A-70400 Includes Set of 4</p>
	<p>Standard Wingset #A-70200 Includes Hip + Chest</p>	<p>Comfort Cover Disposables #A-70260 1 Case of 5 Sets</p>	<p>Supine Top Pads #A-70435 (3" / 7.6 cm) #A-70434 (2" / 5 cm) Includes Set of 4 (Specify Thickness)</p>
Choose Rail Type	<p>Flex Frame OR Table Clamp Set #A-70510-US, #A-70510-DEN, #A-70510-EUR, #A-70510-UK, #A-70510-Jpn Includes 1 Pair, Specify Rail Type</p>	<p>Flex Frame Armboard Rail #A-70500-US, #A-70500-DEN, #A-70500-EUR, #A-70500-UK, #A-70500-Jpn Includes 2, Specify Rail Type</p>	<p>Easy Lock™ Socket #A-40022 (US), #A-40023 (EU), #A-40024 (UK), #A-40030 (Jpn) Includes 2, Specify Rail Type</p>
What's Required	<p>Comfort Mask Disposables #A-70310 1 Case of 5 Required</p>	<p>Comfort Cover Disposables #A-70260 1 Case of 5 Sets Required</p>	
What's Optional	<p>Prone Spinal Leg Pad #A-70441</p>	<p>Short Accessory Rail #A-70501-US, #A-70501-Jpn, #A-70501-DEN, #A-70501-EU, #A-70501-UK Specify Rail Type</p>	<p>Long Accessory Rail #A-70502-US, #A-70502-Jpn, #A-70502-DEN, #A-70502-EU, #A-70502-UK Specify Rail Type</p>
			<p>Contour Breast Positioner #A-70270</p>