

# Metal Rollerboards are Obsolete!

Unique Allen Transfer System... Better for Patients... Easier on Nurses...



Long Board



Short Board



- **Comfortable:** Soft but sturdy board is easy on patients
  - **No Rollers:** Less traumatic transfer for your patients
  - **Lightweight:** Very easy to use
  - **Fowler-Position Transfer:** Long board allows for transfer in seated position
- No Special Training Requirements

## Buy 5 and Receive 10% Off!

Mention Promo G55WS.  
Offer Valid  
March 20 - May 05, 2010.

# Buy 5 and Receive 10% Off!

Mention Promo G55WS.  
Offer Valid  
March 20 - May 05, 2010.

## Allen Patient Transfer Board


G55WS-A-83000 Long Patient Transfer Board  
71.5" x 16.5" (180 cm x 42 cm)

G55WS-A-83001 Short Patient Transfer Board  
40" x 16.5" (100 cm x 42 cm)

## Ultra-Lightweight!

Short Board is 5 lbs. (2.3 kg)

Long Board is 6.3 lbs. (2.8 kg)



## How it Works:

The two low friction surfaces allow the transfer cover to rotate around the transfer board **without metal rollers.**

## Make Patient Transfers Easier for Your Staff.

Visit [www.AllenMedical.com](http://www.AllenMedical.com)  
to see all current promotions.



A Hill-Rom Company

One Post Office Square Acton, MA 01720

PRSR STD  
U.S.POSTAGE  
**PAID**  
BATESVILLE, IN  
PERMIT NO. 7714

RETURN SERVICE REQUESTED

(800) 433-5774 • Fax (978) 263-8846

We accept VISA® and MasterCard®. All of our products are latex-free.

Allen Medical Systems, Inc. reserves the right to make changes without notice in design, specifications, and models. The only warranty Allen Medical Systems, Inc. makes is the express written warranty extended on the sale of its products. Products may be comprised of new, refurbished, and recycled materials. All refurbished and recycled materials meet our stringent new product specifications.

© 2010 Allen Medical Systems Inc. All Rights Reserved  
D-770179-B3 February 26, 2010