


Metal Rollerboards are Obsolete!

Allen Transfer System... Better for Patients... Easier on Nurses...



- **Comfortable:** Soft but sturdy board is easy on patients
- **No Rollers:** Less traumatic transfer for your patients
- **Lightweight:** Very easy to use
- **Fowler-Position Transfer:** Long board allows for transfer in seated position
- **Resilient and Sturdy** - More durable than existing boards on the market*

No Special Training Requirements

Buy 5 and Receive 10% Off!

Mention Promo G71WS.
Offer Valid
Sept. 01 - Oct. 05, 2010.



A Hill-Rom Company

Buy 5 and Receive 10% Off!

Mention Promo G71WS.
Offer Valid
Sept. 01 - Oct. 15, 2010.

Allen Patient Transfer Board

- G71WS-A-83005 Long, Narrow Patient Transfer Board
71.5" x 16.5" (180 cm x 42 cm)
- G71WS-A-83000 Long, Wide Patient Transfer Board
71.5" x 20.5" (180 cm x 52 cm)
- G71WS-A-83001 Short, Narrow Patient Transfer Board
40" x 16.5" (100 cm x 42 cm)



How it Works:

The two low friction surfaces allow the transfer cover to rotate around the transfer board **without metal rollers.**

Ultra-Lightweight!

- Short, Narrow Board is 5 lbs. (2.3 kg)
- Long, Narrow Board is 6.3 lbs. (2.8 kg)
- Long, Wide Board is 7 lbs. (3.2 kg)

Make Patient Transfers Easier for Your Staff.

Visit www.AllenMedical.com
to see all current promotions.



A Hill-Rom Company

One Post Office Square Acton, MA 01720

FIRST CLASS MAIL
U.S. POSTAGE
PAID
BATESVILLE, IN
PERMIT NO. 7714

RETURN SERVICE REQUESTED

(800) 433-5774 • Fax (978) 263-8846

We accept VISA® and MasterCard®. All of our products are latex-free.

Allen Medical Systems, Inc. reserves the right to make changes without notice in design, specifications, and models. The only warranty Allen Medical Systems, Inc. makes is the express written warranty extended on the sale of its products. Products may be comprised of new, refurbished, and recycled materials. All refurbished and recycled materials meet our stringent new product specifications.

*Based on comparison to Samarit Rollbord™

© 2010 Allen Medical Systems Inc. All Rights Reserved
D-770179-B5 August 17, 2010