


Astrus Elite & Astrus Elite Silver Surfaces

Effective Pressure Management Our Astrus Elite™ and Astrus Elite Silver™ Surfaces have a four-layer foam composite construction that protects patients four ways. A top layer of conforming foam molds to your patients' anatomy. The pressure displacement foam middle layer transfers pressure off the bony prominences, and the closed cell foam bottom layer prevents bottoming out. The Anti-Shear™ cover protects patients' skin. The Astrus Elite Silver ticking kills odor-causing bacteria to protect the surface of the pad. Offered with integral gel layer (gray bottom) or without (black bottom.) Available in 3" (7.6 cm), and 4" (10 cm.) Velcro Kit included. Non-conductive.

Astrus Elite Silver offers the same superior pressure management as the Astrus Elite, but now with silver embedded in the cover of the pad!

Astrus Elite & Astrus Elite Silver OR Table Surface: Pressure Mapping Results



Astrus Elite OR Table Surfaces

Astrus Elite - Ordering Info			Astrus Elite Gel - Ordering Info			
3" Pad	4" Pad	5" Pad	3.25" Pad	4.25" Pad	5.25" Pad	
Product No.	Product No.	Product No.	Product No.	Product No.	Product No.	
Berchtold Operon B810, B710	H-303-12-100	H-403-12-100	H-503-12-100	H-313-12-100	H-413-12-100	H-513-12-100
Berchtold Operon D850, D750	H-303-12-110	H-403-12-110	H-503-12-110	H-313-12-110	H-413-12-110	H-513-12-110
Maquet 1130	H-303-10-110	H-403-10-110	H-503-10-110	H-313-10-110	H-413-10-110	H-513-10-110
Maquet 1130 Split Leg	H-303-10-100	H-403-10-100	H-503-10-100	H-313-10-100	H-413-10-100	H-513-10-100
Maquet 1130 Orthopedic	H-303-10-130	H-403-10-130	H-503-10-130	H-313-10-130	H-413-10-130	H-513-10-130
Maquet Alphamaxx	H-303-10-150	H-403-10-150	H-503-10-150	H-313-10-150	H-413-10-150	H-513-10-150
Maquet Alphastar 1132	H-303-10-120	H-403-10-120	H-503-10-120	H-313-10-120	H-413-10-120	H-513-10-120
Midmark 7100	H-303-06-100	H-403-06-100	H-503-06-100	H-313-06-100	H-413-06-100	H-513-06-100
Midmark Chic CST	H-303-06-120	H-403-06-120	H-503-06-120	H-313-06-120	H-413-06-120	H-513-06-120
Midmark Montreal	H-303-06-130	H-403-06-130	H-503-06-130	H-313-06-130	H-413-06-130	H-513-06-130
Getinge / Champagne 4900B, 4900E, 5100B, 5100E	H-303-04-100	H-403-04-100	H-503-04-100	H-313-04-100	H-413-04-100	H-513-04-100
Siemens Universal	H-303-13-100	H-403-13-100	H-503-13-100	H-313-13-100	H-413-13-100	H-513-13-100
Skytron 5001, 6001, 6500, 6600	H-303-05-130	H-403-05-130	H-503-05-130	H-313-05-130	H-413-05-130	H-513-05-130
Skytron 3100	H-303-05-100	H-403-05-100	H-503-05-100	H-313-05-100	H-413-05-100	H-513-05-100
Skytron 3500	H-303-05-120	H-403-05-120	H-503-05-120	H-313-05-120	H-413-05-120	H-513-05-120
Skytron 3600	H-303-05-190	H-403-05-190	H-503-05-190	H-313-05-190	H-413-05-190	H-513-05-190
Skytron 6700, 6701	H-303-05-200	H-403-05-200	H-503-05-200	H-313-05-200	H-413-05-200	H-513-05-200
Steris-Amsco 1080, 2080	H-303-03-100	H-403-03-100	H-503-03-100	H-313-03-100	H-413-03-100	H-513-03-100
Steris-Amsco 3080, 3085	H-303-03-110	H-403-03-110	H-503-03-110	H-313-03-110	H-413-03-110	H-513-03-110
Steris C-Max	H-303-03-150	H-403-03-150	H-503-03-150	H-313-03-150	H-413-03-150	H-513-03-150
Steris Surgi-graphic	H-303-03-170	H-403-03-170	H-503-03-170	H-313-03-170	H-413-03-170	H-513-03-170

Astrus Elite Silver OR Table Surfaces



Astrus Elite Silver - Ordering Info			Astrus Elite Silver Gel - Ordering Info			
3" Pad	4" Pad	5" Pad	3.25" Pad	4.25" Pad	5.25" Pad	
Product No.	Product No.	Product No.	Product No.	Product No.	Product No.	
Berchtold Operon B810, B710	H-30A-12-100	H-40A-12-100	H-50A-12-100	H-31A-12-100	H-41A-12-100	H-51A-12-100
Berchtold Operon D850, D750	H-30A-12-110	H-40A-12-110	H-50A-12-110	H-31A-12-110	H-41A-12-110	H-51A-12-110
Maquet 1130	H-30A-10-110	H-40A-10-110	H-50A-10-110	H-31A-10-110	H-41A-10-110	H-51A-10-110
Maquet 1130 Split Leg	H-30A-10-100	H-40A-10-100	H-50A-10-100	H-31A-10-100	H-41A-10-100	H-51A-10-100
Maquet 1130 Orthopedic	H-30A-10-130	H-40A-10-130	H-50A-10-130	H-31A-10-130	H-41A-10-130	H-51A-10-130
Maquet Alphamaxx	H-30A-10-150	H-40A-10-150	H-50A-10-150	H-31A-10-150	H-41A-10-150	H-51A-10-150
Maquet Alphastar 1132	H-30A-10-120	H-40A-10-120	H-50A-10-120	H-31A-10-120	H-41A-10-120	H-51A-10-120
Midmark 7100	H-30A-06-100	H-40A-06-100	H-50A-06-100	H-31A-06-100	H-41A-06-100	H-51A-06-100
Midmark Chic CST	H-30A-06-120	H-40A-06-120	H-50A-06-120	H-31A-06-120	H-41A-06-120	H-51A-06-120
Midmark Montreal	H-30A-06-130	H-40A-06-130	H-50A-06-130	H-31A-06-130	H-41A-06-130	H-51A-06-130
Getinge / Champagne 4900B, 4900E, 5100B, 5100E	H-30A-04-100	H-40A-04-100	H-50A-04-100	H-31A-04-100	H-41A-04-100	H-51A-04-100
Siemens Universal	H-30A-13-100	H-40A-13-100	H-50A-13-100	H-31A-13-100	H-41A-13-100	H-51A-13-100
Skytron 5001, 6001, 6500, 6600	H-30A-05-130	H-40A-05-130	H-50A-05-130	H-31A-05-130	H-41A-05-130	H-51A-05-130
Skytron 3100	H-30A-05-100	H-40A-05-100	H-50A-05-100	H-31A-05-100	H-41A-05-100	H-51A-05-100
Skytron 3500	H-30A-05-120	H-40A-05-120	H-50A-05-120	H-31A-05-120	H-41A-05-120	H-51A-05-120
Skytron 3600	H-30A-05-190	H-40A-05-190	H-50A-05-190	H-31A-05-190	H-41A-05-190	H-51A-05-190
Skytron 6700, 6701	H-30A-05-200	H-40A-05-200	H-50A-05-200	H-31A-05-200	H-41A-05-200	H-51A-05-200
Steris-Amsco 1080, 2080	H-30A-03-100	H-40A-03-100	H-50A-03-100	H-31A-03-100	H-41A-03-100	H-51A-03-100
Steris-Amsco 3080, 3085	H-30A-03-110	H-40A-03-110	H-50A-03-110	H-31A-03-110	H-41A-03-110	H-51A-03-110
Steris C-Max	H-30A-03-150	H-40A-03-150	H-50A-03-150	H-31A-03-150	H-41A-03-150	H-51A-03-150
Steris Surgi-graphic	H-30A-03-170	H-40A-03-170	H-50A-03-170	H-31A-03-170	H-41A-03-170	H-51A-03-170