

Allen Gel Surfaces

- Advanced viscoelastic formulation provides continuous tissue support that will not bottom out.
- Easy-to-clean flexible membrane cover has shear characteristics similar to that of human skin.
- Will not absorb body fluids or odors.
- Radiolucent and latex-free.
- Repairable if punctured.

Overlay Table Pads



Universal 3-Piece OR Table Pad Set
#A-10034



Head Section Gel Pad
9" x 20" x 3/8" (225 x 500 x 10 mm)
#A-10034H



Torso Section [w/cutout]
39" x 20" x 3/8" (990 x 500 x 10 mm)
#A-10034T



Foot Section Gel Pad
19" x 20" x 3/8" (500 x 500 x 10 mm)
#A-10034F



3/4 Length Table Gel Pad
45 1/4" x 20 1/2" x 3/8" (1150 x 520 x 10 mm)
#A-10060



Full Length Table Gel Pad
71" x 20 1/2" x 3/8" (1800 x 520 x 10 mm)
#A-10061

Head Positioners



Adult Head Gel Donut
8" O.D. x 2" (200 x 45 mm)
#A-10031



Pediatric Head Gel Donut
3 1/2" O.D. x 3/4" (90 x 20 mm)
#A-10103



Adult Gel Horseshoe
8" O.D. x 3" (200 x 75 mm)
#A-10104

Positioners



Small Gel Positioner
12" x 4 3/4" x 2 3/8" (300 x 120 x 60 mm)
#A-10073



Large Gel Positioner
14" x 5" x 2 3/4" (360 x 125 x 70 mm)
#A-10074



Chest Roll Gel Positioner
12" x 4" x 3" (305 x 100 x 75 mm)
#A-10076



Chest Roll Gel Positioner w/ Foam Core
16" x 6" x 6" (405 x 150 x 150 mm)
#A-10097



Chest Roll Gel Positioner
16" x 4" x 3" (405 x 100 x 75 mm)
#A-10078



Chest Roll Gel Positioner w/ Foam Core
20" x 6" x 6" (510 x 150 x 150 mm)
#A-10098



Chest Roll Gel Positioner
20" x 4" x 3" (510 x 100 x 75 mm)
#A-10080

Nerve and Heel Protectors



Small Ulnar Nerve Gel Pad w/ Velcro Strap
15" x 4 1/2" x 5/16" (380 x 115 x 8 mm)
#A-10102



Large Ulnar Nerve Gel Pad w/ Velcro Straps
18" x 8 1/4" x 5/16" (460 x 210 x 8 mm)
#A-10066



Foot-Heel Gel Pad w/ Velcro
4 1/8" x 3 1/2" x 1 3/4" (105 x 90 x 43 mm)
#A-10067



Clamshell Heel Gel Pad w/ Velcro Straps
9" x 7 7/8" x 7/8" (230 x 200 x 20 mm)
#A-10068



Polymer Gel Boot Pads
19 7/8" x 5 3/4" x 3/8" (505 x 145 x 10 mm)
#A-10026

Armboard and Toboggan Pads



Small Armboard Gel Pad
15 1/2" x 4 1/2" x 3/8" (395 x 115 x 10 mm)
#A-10105



Toboggan Gel Pad
15 1/2" x 10" x 3/8" (395 x 245 x 10 mm)
#A-10029



Large Armboard Gel Pad
23 5/8" x 5" x 3/8" (600 x 125 x 10 mm)
#A-10032



Allen Armboard Gel Pad
23 5/8" x 8" x 3/8" (600 x 200 x 10 mm)
#A-10035


Knee Crutch Pads



Small Knee Crutch Gel Pad
12 5/8" x 8 5/8" x 3/8" (320 x 220 x 10 mm)
#A-10070



Medium Knee Crutch Gel Pad
15 1/2" x 10" x 3/8" (395 x 245 x 10 mm)
#A-10029



Hip/Large Knee Crutch Gel Pad
20 1/2" x 20 1/2" x 3/8" (520 x 520 x 10 mm)
#A-10033

Other

Wilson Frame Gel Pad	#A-10101	29" x 10" x 3/8" (735 x 255 x 10 mm)
Repair Kit	#Kit002	4" x 4" (101 x 101 mm)

Protect and Pamper Your Patients with Luxurious Allen Polymer Gel Pads

Note: All OR Table Surfaces, Armboard Surfaces, Stretcher Surfaces, Gel Pads (including Vacu-Gel, Vacu-Form and FlexiForm), and Foam Positioners, returned within thirty (30) days of the invoice date are subject to a 20% restocking fee. This applies to any pads or surfaces returned in new and unused condition.