

GYN/URO/LAPRO



Stirrup Selection & Comparison Chart

	Product	Lift Mechanism	Lithotomy Adjustment (From Sterile Field)	Lithotomy Range (From Sterile Field)	Abduction Adjustment (From Sterile Field)	Abduction Range (From Sterile Field)	Lithotomy Indicator	Patient Weight Limit	Mechanism Cover	Protective Fin Design
Best	Ultrafin® #O-UFAS	Y	Y All Positions in Range	+84, -33	Y All Positions in Range	+25, -9	Y (on rod)	800 lbs. (363 kg)	N	Y
	Yellofin® Elite #O-YFES	Y	Y All Positions in Range	+84, -33	Y All Positions in Range	+25, -9	Y (on rod)	500 lbs. (227 kg)	N	Y
	Yellofin® #O-YFASI	Y	Y All Positions in Range	+84, -33	Y All Positions in Range	+25, -9	Y (on rod)	350 lbs. (159 kg)	N	Y
Better	PAL Pro® #A-10052	Y	Y All Positions in Range	+90, -22	Y All Positions in Range	+25, -9	Y (on handle)	350 lbs. (159 kg)	Y	N
	PAL® #A-10023	Y	Y All Positions in Range	+90, -22	Y Fixed Abduction Track	Fixed Track 2° for 10° lithotomy	Y (on handle)	350 lbs. (159 kg)	Y	N
	Yellofin® (without Lift-Assist) #O-YFSI	N	Y All Positions in Range	+90, -35	Y All Positions in Range	+25, -9	Y (on rod)	350 lbs. (159 kg)	N	Y
Good	SAM III™ #A-10024	N	Y 10 Positions in Range	+90, -26	Y Fixed Abduction Track	Fixed Track 2° for 10° Lithotomy	Y (on rod joint)	350 lbs. (159 kg)	N	N
	Adjustable Length #A-10008-F	N	N	N	N	N	N	350 lbs. (159 kg)	N	N
	Basic Stirrups #A-10001-F	N	N	N	N	N	N	350 lbs. (159 kg)	N	N
Pediatric	Junior PAL #A-10056	N	Y All Positions in Range	+90, -35	Y All Positions in Range	+25, -9	Y (on rod)	250 lbs. (113 kg)	N	N
	Little PAL #A-10057	N	Y All Positions in Range	+90, -35	Y All Positions in Range	+25, -9	Y (on rod)	160 lbs. (73 kg)	N	N

Good
 Better
 Best

*O = Optional

* Extended Warranty Program For Stirrups - U.S. Only

Keep your stirrups functioning like new for 4 full years! Purchase the extended warranty and Allen Medical will extend the original warranty for an additional two years (4 years total) for stirrup hardware. This program ensures that the stirrup functions properly for 4 years from the date of purchase at no additional cost. Warranty Extension continues the original stirrup warranty for two years. This does not include boot pads purchased in conjunction with the stirrup or as part of the original stirrup. See page 113 for more information on Allen Medical System's Warranty Program.