



Product Shown on Top of Table Pad for Clarity

## Lateral Supports

Grip patient's front and back for secure lateral positioning. The supports slide in and out on positioning bar and lock securely in place with a quick twist of the knob. Each soft pressure management padded support has limited rotation to conform to the patient's anatomy.

#O-LBL Lateral Supports (1 Pair)

What's Required



Easy Lock Blade Clamp  
#A-40040 (US), #A-40041 (EU),  
#A-40042 (UK), #A-40043 (Jpn)  
2 Required

What's Optional



Accessory Cart  
#A-30010 (US),  
#A-30010-E (EU),  
#A-30010-UK (UK)  
Specify Rail Type



Replacement Pad  
#O-AMA70X1.3S  
(3" x 9" x 1 1/4")  
(7.6 cm x 23 cm x 3.1 cm)



## Vacu-Gel® & Vacu-Form® Lateral Positioners

Stable and comfortable lateral positioning. Exclusive integral gel layer provides extra protection against pressure ulcer development in Vacu-Gel Positioner. Repair kit included.

- Easy-to-operate valve uses 1/4" (0.6 cm) suction line
- Baffles control bead migration
- Large size, 40" x 30" (102 cm x 76 cm)
- Axillary cutout
- Flap to clip instruments safely

#H-018-00-201 Vacu-Gel Lateral Positioner

#H-008-00-201 Vacu-Form Lateral Positioner

## Legholder Disposable Pads

CE



The Legholder Disposable Pad is compatible with all Allen leg holder systems. These pads utilize two special foam materials which are laminated together. The outer pad is a dense, slightly compressible material which allows for high clamping pressure to stabilize the leg. The inner pad is a soft, fully compressible material for gentle, non-abrasive contact with the skin. The leg pad secures with simple Velcro fasteners. Length is 37 3/4" (95.9 cm) and width is 3 3/8" (8.6 cm).

#R-S4133-V Legholder Pad (Case of 12)



## Deluxe Arthroscopic Legholder System

CE

Smooth, uni-directional locking mechanism enables the caregiver to set rigid fixation with one hand operation. The hinged upper brace accommodates all size patients while the structurally enhanced design enables the staff to apply pressure both medially and laterally. The ratchet-style securement strap easily adjusts for precise amount of compression. Holds legs up to 35" (89 cm) in circumference. Supports up to a 500 lbs. (227 kg) patient.

#A-10015 Deluxe Arthroscopic Legholder System

What's Included	 Deluxe Arthroscopic Legholder System #A-10015	 Leg Holder Disposable Pad #R-S4133-V 3 Included
What's Required	 Easy Lock Socket #A-40022 (US), #A-40023 (EU), #A-40024 (UK), #A-40030 (Jpn) 1 Required, Specify Rail Type	 Leg Holder Disposable Pad #R-S4133-V Case of 12
What's Optional	 Arthroscopic Well-Leg Holder #O-AWLH	 Foam Well-Leg Holder #H-009-00-120
		 Accessory Cart #A-30010 (US), #A-30010-E (EU), #A-30010-UK (UK) Specify Rail Type



## Classic Arthroscopic Legholder System

CE

The Arthroscopic Legholder immobilizes the thigh and allows for lower leg manipulation during knee arthroscopy. A removable slotted post facilitates the placement of the patient's leg into the holder while the stabilizing bar restricts upper leg movement during knee manipulation.

#R-734-33 Classic Arthroscopic Legholder System

What's Included	 Classic Arthroscopic Legholder System #R-734-33	 Leg Holder Disposable Pad #R-S4133-V 3 Included
What's Required	 Easy Lock Socket #A-40022 (US), #A-40023 (EU), #A-40024 (UK), #A-40030 (Jpn) 1 Required, Specify Rail Type	 Leg Holder Disposable Pad #R-S4133-V Case of 12
What's Optional	 Arthroscopic Well-Leg Holder #O-AWLH	 Foam Well-Leg Holder #H-009-00-120
		 Replacement Slide Bar #R-734-34