



Unobstructed Posterior and Lateral Access



Lift-Assist Pistons Make Positioning Patients Easy



Chair Lies Flat for Clear Anesthesia Access

Lift-Assist™ Beach Chair

CE

Make Positioning Easy. Use Lift-Assist One-Handed Positioning.

Integral Clamps!

The Lift-Assist Beach Chair is the latest and most advanced option from Allen Medical Systems for positioning during shoulder procedures. The product offers unobstructed posterior access, and repositioning the patient is quick and easy. The lift-assist design allows for nearly effortless positioning as the pistons support most of the patient's weight.

- **Choice in Positioning** - Positioner attaches next to either break of the surgical table seat section
- **Lift-Assist** - Dramatically reduces the effort to lift heavy patients
- **500 lbs. (227 kg) Weight Capacity** - Helps position larger patients
- **Integrated Clamps with Handles** - Makes it easier to set-up
- **Enhanced Width Extension** - 20⁵/₈" - 24⁷/₈" (52.4 cm - 63.2 cm). Allows chair to mount on a wide array of tables.
- **Head Rest Design** - Allows greater movement to accommodate shorter patients
- **Unique Table Pad Integration** - Prevents gapping between the pad and chair
- **Sliding Headrest** - Reduces stress on neck when raising and lowering the patient
- **Expanded Surgical Access** - Large cut-out for excellent surgical access
- **Lightweight** - Under 30 lbs. (14 kg) allowing for a single-person setup
- **Articulation Range** - 0° - 90°

#A-91000 Lift-Assist Beach Chair

What's Required	What's Included	 Lift-Assist Beach Chair #A-91000 (US), #A-91000E (EU), #A-91000D (Denyer), #A-91000UK (UK) Specify Rail Type	 Reusable Pad Set #A-90020 1 Included	 Reusable Security Strap #A-90003 1 Included		
	Counter-Traction	 Counter-Traction Strap #A-90015 1 Required	— AND —	 Easy Lock Socket #A-40022 (US), #A-40023 (EU), #A-40024 (UK), #A-40030 (Jpn) 1 Required, Specify Rail Type		
	Choose a Headrest	 FlexiForm Head Rest #A-90017 1 Required	— OR —	 FlexiForm Head Rest Pad #A-90019 1 Case of 10	 Standard Head Rest #A-90016 1 Required	
	Choose an Arm Support	 LPS Arm Support #A-90002 1 Required	— OR —	 Arm Platform #A-90007 1 Required	— AND —	 Easy Lock Socket #A-40022 (US), #A-40023 (EU), #A-40024 (UK), #A-40030 (Jpn) 1 Required, Specify Rail Type
	What's Optional	 Beach Chair Cart #A-90100 (US), #A-90100E (EU), #A-90100 (UK) Specify Rail Type	 Allen Intraoperative Arm Positioner #A-92000 (US), #A-92000E (EU), #A-92000D (Denyer), #A-92000UK (UK) Specify Rail Type			