




Unobstructed Access



Effortless Positioning



External Rotation and Abduction



Unlimited Range of Motion

Allen Intraoperative Arm Positioner™



Spend Your Time on the Procedure, Not the Positioning!

Unlimited Surgical Site Access in the Sterile Field







The Allen Intraoperative Arm Positioner delivers effortless arm positioning for shoulder surgery in the Beach Chair position, as well as elbow and clavicle procedures. By indexing at the elbow, the device is able to gain complete positioning control of the shoulder joint. The revolutionary patent pending design allows the entire device to be steam sterilized and provides an unprecedented level of access to the surgical site. The device also has the capability of achieving a semi-locked state, allowing the surgeon to simply move the arm where he wants it without having to adjust the tension.

Ideal for use with all open and arthroscopic upper limb procedures.

Simple: No hoses, gases, power cords or hydraulics

- Effortless, fluid positioning for enhanced site access
- Connects over the drape, in the sterile field
- Compact and portable for easy storage
- Single knob releases all joints
- Index at the elbow for complete control
- The only fully autoclavable arm positioner on the market
- Quick and easy setup
- Unlimited access to the surgical site (no bars or foot pedal obstructions)
- No Foot Pedal / No Hoses / No Hydraulics / No Power Cords
- Lightweight, device weighs only 6 lbs. (2.7 kg). No cart required!
- Can support 500 lb. (227 kg) maximum patient weight

#A-92000 Allen Intraoperative Arm Positioner

What's Included	 Allen Intraoperative Arm Positioner #A-92000 1 Included	 Allen Intraoperative Arm Positioner Clamp #A-92002 1 Included	
What's Required	 Allen Intraoperative Arm Positioner Sterile Disposables #A-92001 Case of 10		
What's Optional	 Lift-Assist Beach Chair #A-91000 (US), #A-91000E (EU), #A-91000D (Denyer), #A-91000UK (UK)	 Table-Powered Beach Chair #A-90000	 Beach Chair Cart #A-90100 (US) #A-90100E (EU) #A-90100 (UK)