



C-Flex Clean Cape Disposables CE

Disposable cover to protect the C-Flex from fluids and other contaminants.

#A-70702 C-Flex Clean Cape Disposables (Case of 10)



C-Prone Head Positioner CE

- Prone head positioner
- Vertical and rotational adjustment
- 500 lbs. (227 kg) patient weight capacity
- Fits Allen Flex Frame and Jackson spine frame
- Fold away mirror
- Eyes free from pressure
- Easy airway access

#A-70300 C-Prone Head Positioner

What's Included



C-Prone Head Positioner
#A-70300



Comfort Mask Disposables
#A-70310
Includes 1 Case of 5

What's Required



Comfort Mask Disposables
#A-70310
1 Case of 5



Supine Top Set CE

Radiolucent sections that convert an open spine frame to a solid imaging surface. 500 lbs. (227 kg) patient weight capacity.

#A-70400 Supine Top Set

What's Included



Supine Top Set
#A-70400
4 Included
(30.5 cm (12") top sections)

What's Required



12" Supine Top Set Pads,
2" (5 cm) Thickness
#A-70434
4 Required

OR



12" Supine Top Set Pads,
3" (7.6 cm) Thickness
#A-70435
4 Required



Flex Frame Armboard Rail CE

Quick connect mount for attaching arm supports to a spine frame. 500 lbs. (227 kg) patient weight capacity.

#A-70500 Armboard Rail

Specify Rail Type

#A-70500 (US) #A-70500-JP (Jpn)
#A-70500-EU (EU) #A-70500-DEN (Denyer)
#A-70500-UK (UK)



Flex Frame Short Accessory Rail CE

Mount retractors anywhere on the carbon fiber frame. Rigid mount for attaching any standard rail accessory to a spine frame. Ideal for retractors and arm supports. 500 lbs. (227 kg) patient weight capacity. Length of rail is 6" (15.2 cm).

#A-70501 Short Accessory Rail

Specify Rail Type

#A-70501 (US) #A-70501-DEN (Denyer)
#A-70501-EU (EU) #A-70501-JP (Jpn)
#A-70501-UK (UK)



Flex Frame Long Accessory Rail CE

Mount retractors anywhere on the carbon fiber frame. Double-length mount for attaching multiple accessories to a spine frame. 500 lbs. (227 kg) patient weight capacity. Length of rail is 17" (43.2 cm).

#A-70502 Long Accessory Rail

Specify Rail Type

#A-70502 (US) #A-70502-DEN (Denyer)
#A-70502-EU (EU) #A-70502-JP (Jpn)
#A-70502-UK (UK)