



Integral Rail Clamp

Table Width Extender 4" (10 cm)

Increases the support surface for larger patients. Clips easily with integral locking mechanism to the rail. Can be used individually or in multiples to create a custom platform. Support surface measures 4" x 20" (10 cm x 51 cm). Supports up to a 400 lbs. (181 kg) patient. Includes side rail for mounting stirrup clamps, armboards, and other accessories. Call Customer Service for replacement pad info.

Table Width Extender 4" (10 cm)

US Rail	EU Rail	UK Rail
#O-TWE1 w/ 1" Pad	#O-TWE1E w/ 2.5 cm Pad	#O-TWE1UK w/ 2.5 cm Pad
#O-TWE2 w/ 2" Pad	#O-TWE2E w/ 5 cm Pad	#O-TWE2UK w/ 5 cm Pad
#O-TWE3 w/ 3" Pad	#O-TWE3E w/ 7.6 cm Pad	#O-TWE3UK w/ 7.6 cm Pad
#O-TWE4 w/ 4" Pad	#O-TWE4E w/ 10 cm Pad	#O-TWE4UK w/ 10 cm Pad



Quick Release Heel Cup

The Twister™ Device

An improved "Candy Cane" design with greater leg clearance. Constructed of stainless steel for smooth operation and improved reliability. Strong rotation lock prevents slippage and permits vertical adjustment from 28" to 41" (71 cm to 104 cm). Supports up to a 300 lbs. (136 kg) patient.

#O-LHPCCA The Twister Device (1 Pair)

What's Required	OR
 Easy Lock Socket #A-40022 (US), #A-40023 (EU), #A-40024 (UK), #A-40030 (Jpn) 2 Required, Specify Rail Type	 Tri-Clamp #A-40016 (US), #A-40016E (EU) 2 Required, Specify Rail Type
 Stirrup Cart #A-30015	 Support Pole Pads #O-TSPP 1 Pair
 Heel Cups #O-HC 1 Pair	 Replacement Straps #P-AMA42X0.1Y-A



Cysto-Lift Knee Crutch

CE

Incorporates our patented stirrup technology allowing for simple and easy positioning of the leg during surgery. No need to access the siderail socket to make leg position changes. One convenient handle controls level of lithotomy and abduction/adduction while maintaining the sterile field. The unique "Lift Assist" neutralizes leg weight allowing even the smallest nurse to easily handle individual legs weighing up to 100 lbs. (45 kg)! Doughy soft pad reduces popliteal pressure. Integral siderail clamp makes attachment to the siderail of any cysto table a snap. Boot type stirrups are strongly recommended for longer procedures. Supports up to a 350 lbs. (159 kg) patient.

#O-CLK Cysto-Lift Knee Crutch (1 Pair)

What's Included	What's Required	OR	What's Optional
 Cysto-Lift Knee Crutch #O-CLK	 Easy Lock Blade Clamp #A-40040 (US), #A-40041 (EU), #A-40042 (UK), #A-40043 (Jpn) 2 Required, Specify Rail Type	 Rail Clamp #O-RC2 (US), #O-RC2XUK (EU/UK), #O-RC2XD (DEN), #O-RCS (Swiss), #O-RC29X32MM (Jpn) 2 Required, Specify Rail Type	 Cysto-Lift Knee Crutch Pad #P-AMA06A0.6S
 Uro Catcher System #O-UC4 (US), #O-UC4E (EU), #O-UC4UK (UK) (See Page 39 for More Info)	 Stirrup Cart #A-30015		



Pads Included

Simplicity™ Knee Crutch

CE

Provides excellent patient comfort and precise surgical positioning for short duration procedures. Doughy soft pad reduces popliteal pressure. One easy to use handle provides vertical adjustment from 16" to 21" (40.5 cm to 53 cm) and 360° horizontal rotation. Supports up to a 400 lbs. (181 kg) patient. Includes Velcro strap for security.

#O-KCP Simplicity Knee Crutch (1 Pair)

Choose A Clamp	OR
 Easy Lock Socket #A-40022 (US), #A-40023 (EU), #A-40024 (UK), #A-40030 (Jpn) 2 Required, Specify Rail Type	 Tri-Clamp #A-40016 (US), #A-40016E (EU) 2 Required, Specify Rail Type
 Replacement Pads #O-KCPP 1 Pair	