



PAL® Stirrups

CE

PAL Stirrups deliver many of the great features of the PAL Pro but with controlled abduction. The abduction safety zone is estimated to be up to 12.5° from midline and the PAL Stirrups abduct at 12.5° as they are elevated.

- Power-assist mechanism is protected by rigid, durable cover
- Indexed rods to aid in duplicating positioning

#A-10023 PAL Stirrups (1 Pair)

What's Required			

Note: Stirrups are Sold in Pairs.



SAM™ III Stirrups

CE

Easy, repeatable lithotomy positioning at an economical price. Just lift and the SAM III Stirrup clicks into place.

- Ten lithotomy levels, indexed rods and auto-abduction provide repeatable leg symmetry
- Slim profile and compact rods do not obstruct work area
- Lightweight design for easy handling

#A-10024 SAM III Stirrups (1 Pair)

What's Required			

Adjustable Length Stirrups

Adjustable rod permits independent adjustment of length and knee flexion. Telescopes from 23.5" to 38" (60 cm to 96.5 cm). Patented Allen adjustable clamping mechanism. Molded boot floats to minimize pressure on the calf.

#A-10008-F Adjustable Length Stirrups (1 Pair)



What's Required		CHOOSE		OR	

Basic Stirrups

CE

Original Allen boot stirrup with patented adjustable clamping assembly locks motion in four directions with single handle for easy, one-step adjustment. Basic Stirrups are economical, and the molded boot floats to minimize pressure on the calf.

#A-10001-F Basic Stirrups (1 Pair)



What's Required		CHOOSE		OR	