

Arm and Hand Tables Selection and Comparison Chart

Allen Medical Systems offers a hand table for every need. All of the tables are radiolucent with outstanding durability, backed by a two-year warranty.


- The Heavy Duty tables are rock steady, perfect for use with microscopes.
- The carbon fiber, Carbon Lights™ tables have no frame or legs on the rectangular and minor procedure table to interfere with radiographic images.
- The Basic tables utilize a phenolic board with a perimeter frame to allow for radiographic imaging.



Rectangular



Hourglass



Minor Procedure

		Heavy Duty	Carbon Lights	Basic
		#R-1034-14A	#A-30450-C	#O-AHTR
Rectangular	Weight	9.3 lbs. (4.2 kg)	5.8 lbs. (2.6 kg)	8.8 lbs. (4 kg)
	Dimensions	18" x 29" (46 cm x 73.7 cm)	15" x 30" (38 cm x 76 cm)	14" x 32" (35.5 cm x 81 cm)
	Material	Carbon Fiber	Carbon Fiber	Phenolic
	Height Adjustable	Yes	Yes	Yes
	Legs	Yes (inner leg only)	Optional	Integral
	Post Type	Square	Round	Blade
	Clamp Type	Included (US, EU and Japan rails only)	Universal Clamp	Rail Clamp
	Pad	1" (2.5 cm) Standard Pad	Sold Separately	2" (5 cm) Standard Pad

		Heavy Duty	Carbon Lights	Basic
		#R-1034-11A	#A-30650-C	#O-AHTH
Hourglass	Weight	11.8 lbs. (5.4 kg)	9.2 lbs. (4.2 kg)	11 lbs. (5 kg)
	Dimensions	23" x 38" (58 cm x 96.5 cm)	16" x 34" (41 cm x 86 cm)	14" x 32" (35.5 cm x 81 cm)
	Material	Carbon Fiber	Carbon Fiber	Phenolic
	Height Adjustable	Yes	Yes	Yes
	Legs	Yes (2)	Integral	Integral
	Post Type	Square	Round	Blade
	Clamp Type	Included (US, EU and Japan rails only)	Integral (US and EU)	Rail Clamp
	Pad	1" (2.5 cm) Standard Pad	Sold Separately	2" (5 cm) Standard Pad

		Carbon Lights™	Basic
		#A-30410-C	#O-AHMPT
Minor Procedure	Weight	3.4 lbs. (1.5 kg)	7.4 lbs. (3.4 kg)
	Dimensions	15" x 24" (38 cm x 61 cm)	14" x 26" (35.5 cm x 66 cm)
	Material	Carbon Fiber	Phenolic
	Height Adjustable	Yes	Yes
	Legs	Legless	Legless
	Post Type	Round	Blade
	Clamp Type	Universal Clamp	Rail Clamp
	Pad	Sold Separately	2" (5 cm) Standard Pad