



Unobstructed Lateral Access



Chair Lies Flat

Table-Powered Beach Chair

CE

Use the Power of Your OR Table.

The Table-Powered Beach Chair easily attaches to the leg section of your surgical table and positions the patient using the power of the table! Supports up to a 400 lbs. (181 kg) patient. Effortlessly positions patients throughout the 0° - 90° articulation range.

- **Power Lift:** Utilizes the power of the surgical table to raise and lower the chair
- **400 lbs. (181 kg) Weight Capacity:** Accommodates most patients
- **Sliding Headrest:** Reduces stress on the neck when raising and lowering the patient
- **Width Extension:** Adjustable width - fits all surgical tables 21 $\frac{3}{8}$ " to 23" (54.3 cm - 58.4 cm) wide
- **Head Rest Design:** Allows greater movement to accommodate shorter patients
- **Unique Table Pad Integration:** Prevents gapping between the pad and chair
- **Sliding Headrest:** Reduces stress on neck when raising and lowering the patient
- **Expanded Surgical Access:** Large cut-out for excellent surgical access
- **Articulation Range:** 0° - 90°

#A-90000 Table-Powered Beach Chair

What's Required	What's Included	 Table-Powered Beach Chair #A-90000	 Reusable Pad Set #A-90020 1 Included	 Reusable Security Strap #A-90003 1 Included		
	What's Required	 Beach Chair Clamp #A-90013 (US), #A-90014 (EU/UK) 2 Required Specify Rail Type				
	Counter Traction	 Counter-Traction Strap #A-90015 1 Required	OR	 Counter-Traction Brace #A-90001 1 Required	AND	 Easy Lock Socket #A-40022 (US), #A-40023 (EU), #A-40024 (UK), #A-40030 (Jpn) 1 Required, Specify Rail Type
	Choose a Headrest	 FlexiForm Head Rest #A-90017 1 Required		 FlexiForm™ Head Rest Pad #A-90019 1 Case of 10	OR	 Standard Head Rest #A-90016 1 Required
	Arm Support	 LPS Arm Support #A-90002 1 Required	AND	 Easy Lock Socket #A-40022 (US), #A-40023 (EU), #A-40024 (UK), #A-40030 (Jpn) 1 Required, Specify Rail Type		
	What's Optional	 Beach Chair Cart #A-90100 (US), #A-90100E (EU), #A-90100 (UK) Specify Rail Type		 Allen Intraoperative Arm Positioner #A-92000 (US), #A-92000E (EU), #A-92000D (Denyer), #A-92000UK (UK) Specify Rail Type		