



Flexible Boot



Clamshell Pad Shown

Ultrafin® Stirrups

CE

Easily Accommodates Patients with Larger Calves!

Note: Stirrups are Sold in Pairs.

The Ultrafin Stirrup is an upgrade of the Yellofin Stirrup. The reinforced rod and joint assembly and the flexible boot make the system capable of positioning a patient weighing up to 800 lbs. (363 kg). Though designed for large patients, the Ultrafin Stirrup is still appropriate for use with the general patient population.

- **800 lbs. (363 kg) Patient Weight Capacity** – Support lithotomy positioning for very large patients.
- **Flexible Boot** – Boot comfortably conforms to the calf shape of the patient while providing support.
 - Fin design helps reduce pressure under the popliteal fossa and superficial peroneal nerve.
 - Flexible boot allows these stirrups to be standardized throughout the facility.
- **Floating Boot** – Self adjusting boot minimizes pressure on the calf, when moving the stirrup.
- **Lift-Assist** – Allows easy movement of the leg when placing it in the desired position.
- **Ultrafin Clamshell Pad** – Completely encapsulates the foot, ankle and calf for comfortable and secure positioning for both large and small patients.
- **Squeeze Grip Handle** – Provides easy intraoperative adjustment without compromising the sterile field. Simply release the handle to secure the leg holder in all directions.
- **Lithotomy Range** – Set the stirrup in any position between +84° to –33° lithotomy.
- **Abduction Range** – Set the stirrup in any position between +25° to –9° abduction.
- **Lithotomy & Length Indicators** – Help ensure precise positioning.

#O-UFAS Ultrafin Stirrups

What's Included	 <p>Ultrafin Stirrups #O-UFAS 1 Pair</p>	 <p>Ultrafin Clamshell Pad #O-AMA65S0.6S 1 Pair</p>
What's Required	 <p>Easy Lock Blade Clamp #A-40040 (US), #A-40041 (EU), #A-40042 (UK), #A-40043 (Jpn) 2 Required, Specify Rail Type</p>	 <p>Safety Drape #A-22000 Case of 10 Pairs</p>
What's Optional	 <p>Stirrup Cart #A-30015</p>	 <p>Uro Catcher System #O-UC4 (US), #O-UC4E (EU), #O-UC4UK (UK) (See Page 39 for More Info)</p>  <p>Extended Warranty #A-XWARUF Available in US Only</p>