



Shoulder Supports CE

Support without direct pressure on the brachial plexus, protecting against nerve damage. Ideal for procedures requiring Trendelenburg. Includes thick foam padding with our flexible cover. Supports up to a 500 lbs. (227 kg) patient.

#A-63500 Shoulder Supports (1 Pair)

What's Required	 Easy Lock Blade Clamp #A-40040 (US), #A-40041 (EU), #A-40042 (UK), #A-40043 (Jpn) 2 Required, Specify Rail Type	OR	 Universal Accessory Clamp #A-40015 (US Only) 2 Required
	<td> What's Optional </td> <td>  Replacement Pad (Left) #63501000 </td> <td>  Replacement Pad (Right) #63502000 </td>		What's Optional



Nissen Strap/Thigh Retractor CE

Secure patients from sliding during steep reverse Trendelenburg procedures such as Nissen Fundoplication. Conductive straps encircle each upper thigh with a broad soft pad. Airplane style buckle permits convenient length adjustment. Pad is 4.5" (11.4 cm) wide and 20" (51 cm) long. Straps are adjustable up to 48" (122 cm).

#O-NTS Nissen Strap (US, 1 Pair)

Specify Rail Type

O-NTS (US)	O-NTSJ (Jpn)
O-NTSE (EU)	O-NTSDEN (Denyer)
O-NTSUK (UK)	

Integral Clamps



Uro Catcher® System CE

Provides excellent fluid containment for hysteroscopy or cysto procedures on any table. Flexible drape support frame provides unobstructed surgeon access during manipulation and springs back to its original state to facilitate optimal fluid drainage. Mounts easily with integral rail clamps to side rails and maintains its funnel shape.

#O-UC4 Uro Catcher System (US)

Specify Rail Type

O-UC4 (US)	O-UC4J (Jpn)
O-UC4E (EU)	O-UC4DEN (Denyer)
O-UC4UK (UK)	

Disposables Sold Separately



Uro Catcher Drapes CE

Open mesh filter allows for fast drainage of fluids. The flexible drain hose is compatible with all suction systems for accurate fluid measurement and metering. Velcro securement eliminates the need to use tape. Durable, latex-free drape has flexible hose. Length of tubing is 64" (162.5 cm) extended, 20" (50.8 cm) compressed. Tubing diameter is 7/8" (2.2 cm).

#O-UC31 Sterile Drape (Case of 10)

#O-UC32 Non-sterile Drape (Case of 10)

Open Mesh



Brachytherapy Drapes for Uro Catcher CE

Our ultra-fine mesh trap is designed to catch even very small particles while allowing fluid to drain through. The flexible drain hose is compatible with all suction systems for accurate fluid measurement and metering. Velcro securement eliminates the need to use tape. Durable, latex-free drape has flexible hose. Length of tubing is 64" (162.5 cm) extended, 20" (50.8 cm) compressed. Tubing diameter is 7/8" (2.2 cm).

#A-50008 Sterile Drape (Case of 10)

#A-50007 Non-sterile Drape (Case of 10)

Ultra-Fine Mesh