



**Scoliosis**



**Prone and Lateral Positioning**



**Pediatric Patients**



**Wingset Disposables**

## Wingset Body Positioners

CE

### Supports That Adjust To Fit The Patient.

The wingsets are prone patient supports that have been optimized for pressure management and adjustability. These supports accommodate varying patient anatomies.

- Allows abdomen to hang free
- Avoids compression of brachial plexus
- Offers sternal decompression
- Adjustable for barrel-chested, pediatric, and scoliosis patients
- Used with low shear, moisture-absorbant comfort cover disposables
- Multi-layer pad design for optimal pressure management
- 500 lbs. (227 kg) patient weight capacity
- Accommodates pediatric to adult patients
- Pads adjust to widths of 10" - 18" (225 mm - 455 mm)
- Fits the Allen Flex Frame
- Fits the Jackson Spine Table

Optional Hip Pad Set with cut-out is ideal for larger patients. The Cut-Out Hip Pad must be used with the Standard Chest Pad Set.

### #A-70200 Standard Wingset System

<b>What's Included</b>	 <p>Standard Wingset #A-70200 Includes Hip + Chest Set</p>	 <p>Comfort Cover Disposables #A-70260 Case of 5 Sets</p>
<b>What's Required</b>	 <p>Comfort Cover Disposables #A-70260 Case of 5 Sets</p>	
<b>What's Optional</b>	 <p>Standard Chest Wingset #A-70203</p>	 <p>Cut-Out Hip Wingset #A-70201</p>