

Minor Procedure

Carbon Lights™ Minor Procedures Table CE



This small and ultra lightweight table works well for procedures such as carpal tunnel releases and AV fistulas. Compatible with any OR table. Has two round 5/8" (1.6 cm) posts. Call Customer Service to inquire about clamps to fit alternate rail sizes.

Note: Pad sold separately.

Weight: Approx. 3.4 lbs. (1.5 kg)
Table Dimensions: 24" x 15" (61 cm x 38 cm)
Patient Weight Capacity: 400 lbs. (181 kg)

Carbon Fiber

#A-30410-C Carbon Lights Minor Procedures Table

#A-30410-CL Carbon Lights Minor Procedures Table w/ Leg

What's Required		CHOOSE		OR	
What's Optional					

Minor Procedure

Basic Carpa™ Table CE



The Carpa Table is a legless design with a rigid self-bracing frame. It has a 14" x 26" (35.5 x 66 cm) radiolucent top for maximum access and weighs a mere 7.4 lbs. (3.4 kg)! The perfect choice for carpal tunnel procedures. Easy drop-in rail attachment. Attaches to any OR table. Includes 1.5" (3.8 cm) pad, for use with 2" (5 cm) OR table pads.

Weight: 7.4 lbs. (3.4 kg)
Table Dimensions: 14" x 26" (35.5 cm x 66 cm)
Patient Weight Capacity: 400 lbs. (181 kg)

Pad Included!

#O-AHMPT Basic Carpa Table

What's Included		
What's Required	Easy Lock Blade Clamp #A-40040 (US), #A-40041 (EU), #A-40042 (UK), #A-40043 (Jpn) 2 Required	
What's Optional	Replacement 3" (7.5 cm) Pad #O-AHMPTP3	


Outer rail requires US rail clamp for accessories

Wrist Fracture System




The Wrist Fracture System provides superior traction without weights. An automatic ratcheting design can be adjusted incrementally to provide up to 21 lbs. (9.5 kg) of tension with a manual override setting for tension relaxation when needed. A unique radiolucent hand plate can be fully rotated in either direction (See figures A and B). The Wrist Fracture System is fully autoclavable.

#R-734-51 Wrist Fracture System



What's Included		
What's Required	Simple Clamp #A-40018 (US), #40018-UK (UK), #A-40018-DEN (Denyer) 1 Required	Heavy Duty Hourglass Table #R-1034-11A 1 Required
What's Optional	Sterile Mesh Finger Trap Disposables See Page 71	The Equalizer #R-734-05

The Equalizer CE



Used when three or more digits are placed in traction, The Equalizer distributes the weight of the arm equally between each digit. This stainless steel chain and pulley system is fully autoclavable and comes with an attaching clip. The Equalizer can also be used with the Wrist Fracture System and Weightless Shoulder Suspension System.

#R-734-05 The Equalizer

What's Required	Sterile Mesh Finger Trap Disposables See Page 71	
What's Optional	Easy Irrigation Tower #O-LPA	Wrist Arthroscopy System #R-734-10