

Allen® Traction Boot

CE

Active Heel™ Technology Provides Superb Traction to Prevent Foot Slip

Call Customer Service for Pricing

The Allen Traction Boot with Active Heel™ Technology provides superb traction results during orthopaedic procedures. This patented technology can be used with a wide range of foot sizes – from size 5 female to size 14 male – to securely hold the patient anatomy and sustain traction during fracture, trauma and arthroscopic procedures. The unique design accommodates either the left or the right foot and provides integral handles for quick and easy patient positioning. An intuitive ratchet mechanism helps provide precise traction. The quick release mechanism provides quick release of traction and enables the patient anatomy to be easily removed from the boot. The boot disposable helps maintain skin integrity, patient comfort and cleanliness.


- Fits size range 5 (female) to 14 (male)
- One boot fits left or right foot
- Active Heel Technology
- Handles for easy foot manipulation



Integral Handles Provide Easy Patient Positioning

#A-93000 Allen Traction Boot

#A-93001 Traction Boot Disposable (not shown)



Unique Design Accommodates Left or Right Foot



Ratcheting Mechanism Provides Precision



Active Heel™ Technology Accommodates Many Foot Sizes



Easy to Use Buckles