

Allen Intraoperative Arm Positioner

The Allen Intraoperative Arm Positioner is the smallest, lightweight arm positioner on the market – for any shoulder procedure in the Beach Chair position. It attaches over the drape and provides effortless positioning of the arm intraoperatively - at a fraction of the price of the Spider™ and the McConnell™ arm positioners. Additionally, the disposable is more intuitive and competitively priced.

Specifications

- Fits Accessory Rail Types: US, EU, UK, DEN, JP
- Device weighs 6 lbs. (2.7 kg)
- 500 lbs. (227 kg) maximum patient weight



Intraoperative Arm Positioner
(A-92000)



Arm Positioner Clamp
(A-92002)




Intraoperative Arm Positioner -
Sterile Disposables (Case of 10)
(A-92001)

Optional Accessories

The Allen Intraoperative Arm Positioner works well with the Allen Lift-Assist Beach Chair or the Allen Table-Powered Beach Chair.



Allen Lift-Assist
Beach Chair (A-91000)



Allen Table-Powered
Beach Chair (A-90000)



Allen Universal
Head Positioner (A-90022)

Allen Medical Systems, A Hill-Rom Company
One Post Office Square, Acton, MA 01720

www.AllenMedical.com

United States 800-433-5774 / 978-266-4200
International +44 (0) 870 112 0769

Allen is a registered trademark of Allen Medical Systems, Inc.
Spider is a trademark of TENET Medical Engineering, Inc.
McConnell is a trademark of McConnell Orthopedic Manufacturing Company.

© 2010 Allen Medical Systems, A Hill-Rom Company. All Rights Reserved.
D-770395-A5 October 25, 2010

Allen®

Allen®



Allen® Intraoperative Arm Positioner™


Unlimited Surgical Site Access In The Sterile Field

ALL PRODUCTS ARE CLASS 1 MEDICAL DEVICES

Spend Your Time On The Procedure, Not The Positioning...



Unlimited surgical site access in the sterile field.



Lightweight / Portable



- Connects OVER the drape in the sterile field
- Simple: No hoses, gases, or hydraulics
- Effortless, fluid positioning for enhanced access
- Compact for easy storage
- Lightweight and portable
- Fully autoclavable

Effortless Positioning: Complete range of motion for all shoulder procedures.



Disposable is intuitive and easy to use