

Allen® Easy Catcher™ Disposable System

The All-in-One Fluid Containment System

Available
Spring 2011!

FULLY DISPOSABLE

- No separate band system needed.
- Nothing to lose or misplace between cases.
- No parts to clean or store!



EASY SET-UP / REMOVAL

- Easy Catcher comes fully assembled, and ready to go. Just remove it from the pouch and attach the clamps to the table rails.
- Single piece system: Integral clamps and integrated band
- No disassembly required: The entire system is removed from the table rails, and easily thrown away.



Same benefits as the Allen® Uro Catcher® System:

- Flexible band support frame allows for unobstructed surgical site access. Band snaps back to its original state for optimal fluid drainage.
- Mounts onto OR table side rails.
- Same bag components including: latex-free drape, mesh filter, patient buttocks flap, and collapsible hose.

Allen[®]
A Hill-Rom Company